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MRF technology

Our latest establishment is a production line for pre-fine grinding, polishing and MRF polishing technology. The "Magneto-Rheological Finishing" process ensures a zonal correction of smallest deviations from the nominal radius. This technology allows the production of aspheres as well as surfaces with an accuracy of less than lambda/10 on controlled conditions.

Optical coating

Lens surfaces transmit approx. 96% of the light (due to reflection). For this reason the surface is coated with a thin dielectric film. This coating consists of 5 to 11 layers (in special cases up to 50 layers) and prevents losses on the glass/air (or air/glass) surface. To reduce reflection up to 0.05% the layer has to be adapted to the required wavelength and glass type. Our standard range covers anti-reflective coatings for a wavelength range from 193 nm to 2000 nm.

Quality control

After all production steps the surface quality of our lenses is controlled with a magnifying glass or a microscope. In order to assure the high quality of our products, Sill Optics has been certified according to DIN EN ISO 9001:2015. This motivates all our employees to maintain maximum quality awareness.

Housings

Sill Optics runs its own shop floor for turning and grinding of housings for prototypes and small quantities. We have installed precision turning machines and a 3D-CNC grinding center. Favorable material used is aluminum (RoHS conform), but other materials like titanium are possible as well. After machining all parts receive a black anodized finish or some other varnish to avoid internal reflections.

Assembly

Optical systems (objectives) consist of several lenses that have to be assembled into an objective mount. The lenses must be thoroughly cleaned to avoid any dirt or dust. Attention must be paid to make sure they get centrically set into their mount. Modular systems ensure cost-efficient assembly as well as effective stock keeping. Here, we assemble up to 100,000 objectives per year.

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Additional expertise and custom specific development

Lenses for laser process observation

Observing lenses

Typically laser process observation in done through scan head and scanning lens themselves (1). If there is a difference between observation and working wavelength, a lateral and longitudinal focus shift, caused by monochromatically corrected laser lenses, arises. The scanner aperture could limit the maximum usable aperture of the imaging lens which decreases the resolution.

Alternatively a transparent working plane can be observed from the opposite side (2).

For applications with high resolution specifications and an opaque working plane the observation can be installed next to the f-theta lens (3). A tilt between the optical axis and the orientation of the object plane results. Depending on the depth of field of the vision lens there is only one stripe in focus which can be imagined with high resolution. This lack of definition can be compensated by a suitable adaptor which is mounted between camera and observing lens.

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Optical coherence tomography

Conventional observation lenses do not suffice for three dimensional scan field measuring. Optical coherence tomography is a contactless method and the solution for this kind of process observation.

The basic idea is the superposition of the wavelength of two beam paths. The generated interferogram helps to calculate a distance difference between reference and measuring beam path. If reference distance and distance difference are known the distance of the measuring beam path can be calculated. Other measurements, which are displaced in x- and y direction create a three dimensional map of the scan field. Accuracies in the range of some micrometers can be reached.

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F-theta lenses for polygon scanners

Lenses for mask imaging

and distortion free features.

maskless or direct write and mask projection.

Standard galvanometer scanners are often not fast enough for lasers with extremely high repetition rates. Ultra fast polygon scanners are a highly suited alternative in this case. They are line scanning systems and much faster than galvanometer scanners, which highly decreases the laser processing time. Many existing f-theta lenses are suitable without significant impact on the performance. However, Sill Optics offers custom lenses to exploit the full advantage of the possibilities of polygon scanners.

There are two methods usually used for micro-structuring with lasers,

Flexibility, ease of use and cost effectiveness are the features of the

direct laser process with DPSS laser for feature sizes on a micrometer

Mask projection systems usually use excimer lasers where a mask

pattern is de magnified by a certain factor to reach feature sizes on a

sub-micrometer scale. Additional benefits are good depth uniformity

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Sill Optics designs telecentric lenses for scanning and distortion free imaging on customers' requirements.

Diff erent housing material for UV lenses

Long term exposure of black anodized aluminium housing parts to ultraviolett light can lead to bleaching. Eventually released debris might contaminate the lens surfaces and result in a decreased lifetime on the optical components. Therefore it is possible to buy UV lenses from Sill Optics with high resistant stainless steel housings. Another cost effi cient alternative are colorless anodized aluminum housings with the same price as the standard version.

Dynamic z shifting optics

Besides the well-known f-theta lenses for galvanometer scanners where the scanner is prior to the scan lens, lens systems with adjustable focal length could be used for processing fl at fi elds. These lens systems incorporate a moving lens and a focusing system. The position of the moving lens in respect to the focusing lens has to be synchronized with the scanner movement to avoid a spherical scan fi eld. So it is possible to create an even scan fi eld, but the beam spot diameter depends on its position and changes more over the fi eld in respect to f-theta lenses. Nevertheless there are some applications which utilize the big advantage of working in three dimensions with z shifting optics.

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Alignment turning

Alignment turning is a high precise production technique which enables a minimum tilt between lens axis and optical axis. Production accuracy is about one micrometer with this tool. Before starting the process it is necessary to fix the lens into its housing by curling, sticking or using knurled rings. After that the lens becomes positioned inside the machine so that optical and mechanical axis are exactly in line with each other. With the help of ceramic tools, the outer surfaces as well as the front and rear bearing surfaces are centered. Outside diameter tolerances of just a few microns are producible. If there are high specifications on lens centering, Sill Optics uses this technique for producing modern high precision optics.

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Beam manipulation

Due to the constantly growing requirements for throughput and precision, in many industrial applications beam splitters and beam shapers are in use.

Beam splitters are optical components that split a beam into two or more beams with a certain angular spacing. To realize this, so-called DOEs, diffractive optical elements are often used. DOEs can generate 1D or 2D patterns. The angles between neighboring orders depend on the order numbers and are not linear. This nonlinearity results in non-equidistant spacing of focal points if standard lenses are used. To compensate this scan lenses can be designed in a special way.

Beam shaper transform a single mode Gaussian beam into a beam with uniform energy distribution. Diffractive or aspherical optical elements are used to achieve this. It should be noted that in subsequent optical systems such as beam expanders or f-Theta objectives, more than twice the beam diameter is often required as a free aperture. In addition, the imaging power must be diffraction limited to obtain the beam shape.

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MultiSpot MAShES NextGen3DBat Adalam InScanFoLa SolderScan Sill Optics produces different rotationally symmetric lenses with the philosophy "100% made in Germany". The product portfolio reaches from spherical lenses, aspheres to segmented lenses, domes and rectangular lenses. We are looking forward to your request on a custom specific lens which is designed especially for your application.

Research projects

Sill Optics will always be open for a partnership in research projects. In the past we concluded different projects. Central themes were laser based micro material processing with USP lasers or multispectral lenses for laser processing. Other projects focused topics like e.g. water beams as light guides for laser radiation or compensation elements for the thermal focus shift.

You can find more information about our research projects on our homepage www.silloptics.de/en at company \rightarrow research activities. If you are looking for an industrial partner in optics production for your research project, please feel free to request at Sill Optics.

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General explanations

Part number

The part number is unique for every lens design. Special requests with changes to catalog lenses are assigned a new part number, so repeated purchase with the same specifications is easily possible. The complete part number consists of a 5-digit lens type declaration (e.g. S4LFT = f-theta lens), a 4-digit design identifier number and a 3-digit wavelength and coating specification. Sill reserves the right to make constructional changes in the course of product improvement.

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Datasheets, technical drawings and CAD-files

Individual datasheets are available for each lens via download from our homepage www.silloptics.de . Just go to "Products" on the top menu. Then click on the category "Laser optics" and then the subcategory of the needed lens, e.g. "f-theta lenses". You see a basic explanation about f-theta lenses and on the top a button called "product overview", which leads to an interactive table with all catalog products of the specific category. On the right of each entry, you can find the download-button for the datasheet. Technical drawings, so-called outlines, and 3D CAD-files also available for download.

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Anti-reflective and low-absorption coatings

Our anti-reflective coatings are optimized for a certain wavelength or wavelength ranges. They allow a high transmittance of the laser light and less absorption of energy in the lens for specific wavelengths. Low-absorption coatings are recommended for lasers with a high average power, as they minimize thermal effects. These coatings are only available for fused silica lenses. Beside our standard coatings we also offer customized coatings. For detailed information about damage thresholds, an extra chapter follows on the next pages.

ext.	type	wavelength	specification
/008	anti-reflective	1500 nm - 1600 nm	R<0.25%
/065	broadband	400 nm - 900 nm	R < 0.5 %, avg.
/075	anti-reflective	355 nm	R<0.2%
/081	anti-reflective	1064 nm 532 nm	R<0.2% R<0.25%
/094	anti-reflective	800 nm - 980 nm	R<0.2%
/121	anti-reflective	532 nm	R<0.2%
/126	anti-reflective	1064 nm	R<0.2%
/159	anti-reflective	1850 nm - 1980 nm	R<0.25%
/173	anti-reflective	400 nm - 410 nm	R<0.2%
/199	anti-reflective	255 nm - 266 nm	R<0.2%
/292	low-absorption	515 nm - 545 nm	R<0.2%
/328	low-absorption	1030 nm - 1090 nm	R<0.2%
/373	anti-reflective	420 nm - 480 nm	R<0.2%
/449	low-absorption	900 nm - 1070 nm	R<0.25%
/450	anti-reflective	1000 nm - 1100 nm	R<0.25%
/574	low-absorption	343 nm - 355 nm	R<0.2%

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Design wavelength

All optical systems, especially the coating of the lenses, are designed for a special wavelength or wavelength range. The specifications are demonstrated for a single wavelength and can deviate for a different wavelength in the wavelength transmission band. For the S4LFT4010/292, the specifications are declared for 532 nm but it is possible to use this f-theta objective in a wavelength range from 515 nm to 545 nm without concern.

Absorption and thermal focus shift

Basic problem

With the increasing powers of available lasers, commonly used optical glasses have been pushed to their limits in respect to acceptable thermal effects: Via the exposure of optical glasses to laser radiation, parts of the beam energy are absorbed into the material. This leads to a heating effect with two mayor influences onto optical properties. First, the heating changes the index of refraction of the glass and second, thermal expansion leads to changes in the surface curvatures and therefore changing the refraction of the laser beam. In the application, starting from average laser powers of about 50 W at 1064 nm, a resulting shift of the focal position can lead to decreasing process quality and make online adjustments necessary.

Possible solution

Another option lies in the use of fused silica as lens element material. It is a very resistive glass type which has also a very low absorption coefficient compared to optical glasses. Therefore it is commonly used to minimize thermal effects. Sill also uses special low-absorption coatings to minimize thermal effects further and increase typical damage thresholds.

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Other influences on absorption

Cleanness of the optical components plays also an important role for absorption and therefore thermal effects. Of course there are any larger dust particles (e.g. finger prints) are an extreme absorber of laser radiation. But also time brings micro particles or other contaminations onto lens surfaces. These are not visible to the human eye, but they can be measured.

Measurements show that the absorption increases with the duration to the last cleaning. A new lens cleaning resets the absorption values to the original state.

Therefore regular cleaning of surfaces exposed to the environment is recommended. This effect could not be shown for internal lens elements, thus cleaning them is not necessary!

Please visit our website for further details on cleanness of lenses and lens elements:

www.silloptics.de \rightarrow products \rightarrow Sill Encyclopedia \rightarrow Correct Lens handling and cleaning

Working distance

The working distance is defined as the distance from the front edge of the housing of the lens to the focal or scanning plane of the objective. Be careful to not mix this up with the effective focal length (EFL) of an objective. This is measured from the principle plane, which is a hypothetical plane were the refraction of the complete lens system can be assumed to occur, to the focal plane of the optical system.

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Diffraction value M²

The ability of focusing laser light is defined by ISO standard 11146 and is described by the diffraction value M^2 . This parameter is defined as the ratio of the divergence angle of the laser beam as compared to the divergence angle of an ideal Gaussian beam. An ideal Gaussian beam would provide the smallest possible focus diameter and would have an M^2 value of 1. Sometimes the quality of the laser beam is also described by a parameter K which is the reciprocal of M^2 . The quality of a fiber laser is often defined by the Beam Parameter Product (BPP). This value is given by the product of the diffraction value M^2 and the wavelength λ divided by π . Sill assumes an M^2 values of 1 in all statements about spot sizes. Multiply the spot diameter by the actual laser M^2 value to obtain an actual spot diameter.

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Reversed Mode

In some applications it is necessary to reduce the beam diameter instead of expanding it. Because of the same behavior of light in forward and reversal mode beam expanders with a magnification of x can be turned around to reach a magnification of 1/x. There are also some applications for the reversal mode of f-theta lenses. Due to the changed lens order, a beam expander or f-theta lens, which is normally free of critical back reflections, is not necessarily ghost free in reversed operation.

In the beam expander datasheet there is an information if the lens is free from internal ghosts and suitable for USP lasers for the reversed mode ("no internal ghosts, reversed usage"). If you are not sure if a lens is suitable for your application (forward or reversed usage) feel free to ask. Based on your system data (wavelength, input beam diameter, pulse duration, pulse energy and cw power) it is possible to check if lens and laser match.

Beam diameter d, (1/e²)

The beam diameter is referenced to the drop of beam intensity to the $1/e^2$ point, i.e. to a point where the intensity has fallen to 13.5% of the peak intensity. If the beam diameter is cut at that point, the 13.5% of the total intensity would be lost. If the clipping occurs at the 1.5-times beam diameter, the loss is only approximately 1%.

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Apodization factor

The beam shape and focused spot size after transmission through a lens is strongly dependent on its entrance profile in comparison to the entrance pupil. A common description is given by the truncation ratio T, which is the entrance beam diameter d_L divided by the clear aperture d_{EP} Typical examples are shown in the following sketch: Below a ratio of 0.5, the beam is approximately untruncated. If the entrance beam diameter at $1/e^2$ is equal to the clear aperture of the lens then T=1. Typical applications are located in between those values, as a compromise between low intensity losses or small spot sizes and high costs due to large diameter lenses.

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$$d_{\rm F} = \frac{\lambda \cdot f' \cdot \rm{APO} \cdot M^2}{d_{\rm L}}$$

Spot size (1/e²)

The minimal attainable focal spot size is calculated by the wavelength of the laser multiplied with the focal length of the scan lens, the APO factor and the diffraction parameter M² of the laser divided by the $1/e^2$ beam diameter d_1 .

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beam with $d_1 = 6.00 \text{ mm}$ and $d_1 = 10.0 \text{ mm}$. We assume the use of a f-theta lens S4LFT4010/292 with a frequency doubled Nd:YAG laser at 532 nm and a diffraction value $M^2 = 1.2$. The lens has an effective focal length of f' = 100 mm. Another very important value to determine in addition to the truncation ratio T is the clear aperture or entrance pupil. This is not the clear aperture of the f-theta lens (Ø35mm), but typically the limiting factor is the beam entrance diameter or aperture of the scan system. Assume a very common value of $d_{ep} = 10.0$ mm in this case.

 $T = \frac{6.00 \,\mathrm{mm}}{1.33} = 0.600 \xrightarrow{\mathrm{graph}} \mathrm{APO} \approx 1.33$

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Rayleigh length z

10.0 mm

In optics and especially laser science, the Rayleigh length or Rayleigh range is the distance along the propagation direction of a beam from the waist to the position where the area of the cross section is doubled.

 $d_{\rm F} = \frac{\lambda \cdot f' \cdot \text{APO} \cdot M^2}{d_{\rm L}} = \frac{532 \,\text{nm} \cdot 100 \,\text{mm} \cdot 1.83 \cdot 1.20}{10.0 \,\text{mm}} \approx 11.7 \,\mu\text{m}$

The Rayleigh length is calculated by the focus area multiplied by a factor (depending on the APO factor) divided by the wavelength and the diffraction value M² of the laser.

$z_{\rm R} = \pi \cdot \left(\frac{d_{\rm F}}{2}\right)^2 \cdot \frac{\left(\text{APO}/1.27\right)^2}{2 - M^2}$

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The depth of focus of the scan lens can be estimated by a doubled Rayleigh length. Be aware, that this is just a rough estimation and in many modern applications this value can be too large to still fulfill needed spot diameter requirements.

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Fiber imaging

A large focal spot is an advantage in many applications, like welding. Here, fiber-guided lasers are often used which, in contrast to fiber lasers, have a rather poor beam quality (M² value). As collimator aspheres are used, where Sill offers standard and custom solutions. In cutting heads the focusing is realized by aspheres, in scanning systems an f-theta optic is used.

Using the ratio of focusing focal length f_2 to collimator focal length f_1 provides an easy way to calculate the approximate spot size by the ratio $M = f_2/f_1$. Usually, the fiber core is magnified. In a few applications a demagnification could be used but is limited due to physics.

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$$d_{\rm focus} = M_{\rm NA} \cdot d_{\rm fiber} \approx 587\,\mu{\rm m}$$

 NA_{focus}

The stated calculation is only a rough estimation, that is valid when both focal length are close to each other. A more precise value can be estimated using the numerical apertures:

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Laser induced damage threshold (LIDT)

Principal description

With the increase of more and more powerful laser systems, conventional optics have to withstand these increasing powers. To specify this behavior, the LIDT (Laser Induced Damage Threshold) was induced. It defines a critical laser power or peak fluence, which the lens bulk materials or the coatings cannot withstand and therefore the laser radiation causes permanent changes in the optical characteristics, permanently damaging the lens. With the broad variety of lasers and different thresholds for the LIDT, it is often a very confusing and vague topic. Therefore we will provide some overview and explanations in the following pages. The basis for this information is the international ISO 21254.

Energy density / Fluence

The most common LIDT declaration is given as the max. energy density - also called fluence – which the used material can withstand. The fluence is a unit of laser pulse energy per focal spot area and mostly used for pulsed lasers. For an expressive declaration of the LIDT, not only the value itself, but also the measurement conditions should be given. This includes laser pulse duration, pulse repetition rate and the wavelength.

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When considering cw (continuous-wave) lasers, a definition of the LIDT by a pulse energy does not make any sense. Therefore it is declared as an irradiance, which is described by the laser peak power per effective area (units W/cm² or W/mm²). It is also possible to define the damage threshold by a linear power density (units W/cm or W/mm). If this value is given for the laser, a measurement of the beam size is relevant.

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Although the power is defined as laser pulse energy per pulse duration, a conversion between a LIDT in fluence to a LIDT in intensity is not automatically valid. The reason lies in the different destruction mechanisms of the coatings from various laser types and test conditions.

Dependencies of LIDT

Laser induced damage is dependent on very many factors and therefore it is not possible to guarantee the performance under all possible circumstances. The given conversions in this chapter are therefore just estimates. To ensure a stable and permanent high quality performance of the lenses, we recommend a usage below 50% of the LIDT. Close operation to the LIDT could cause long-term damage. Another important influence originates from tidiness and periodic cleaning, which will be addressed in an extra chapter.

LIDT values depend – beside the pulse energy or peak power - on the laser wavelength λ [nm], pulse repetition rate R[typ.Hz–MHz], pulse duration τ [fs–ms], laser spot diameter $d_{L}[\mu m - cm]$ / laser spot area $A[\mu m^{2} - cm^{2}]$ and the pulse shape. Of course, most of these values only affect the fluency LIDT, as they are pulsed laser values. For cw lasers, only wavelength, pulse shape and sport size are important.

In the following explanations of LIDT dependencies, examples given assume an LIDT of 1 J/cm², measured at 1064 nm with a pulse duration of 1 ns and a repetition rate of 50 Hz.

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Wavelength λ

For any change in the laser wavelength, the LIDT value has to be adjusted. The energy of a photon depends on the wavelength $E = \hbar \omega = 2 \pi c / \lambda$, so shorter-wavelength laser light carries photons with higher energy. This relation also reflects the behavior of the LIDT: Half the wavelength results in the half LIDT value. Great care must be taken with shorter wavelengths. These photons carry much more energy and therefore the actual LIDT value can be much lower than estimated! Thus declarations of the LIDT are meaningfully specified in the wavelength region of their use.

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Pulse duration τ

The dependency on the pulse duration scales with the square root. For example, a 100 ns pulse has an approx. LIDT $\sqrt{100} = 10$ times higher LIDT than that for a 1 ns pulse with comparable pulse energy. This can be counter intuitive! This confusion comes from the separation of the time factor in comparison to the area factor. When calculating the fluency of the laser beam, one does divide the laser pulse energy by the beam area. But at no point in this calculation there is a consideration that laser pulses compress the laser pulse energy in their pulse duration. This means, that a pulse with a shorter pulse duration, but the same overall laser pulse energy, has to compress the same amount of energy into a shorter time interval. This leads to higher peak intensities and therefore a much higher probability of damage. But against intuition, this is not included in the "fluency" but the LIDT value is only given for a specific pulse duration.

In conclusion a shorter laser pulse does decrease the LIDT value as it has higher peak intensities and therefore a higher chance to damage the material. Thus a conversion of the LIDT value is necessary for differing pulse lengths, as it is not directly considered.

Please be aware of this conversion being only valid in the range from around 30 ps to many ms pulse duration. Below the mentioned 30 ps, multi-photon absorption starts to occur, introducing new destruction mechanisms and therefore any optic has lower thresholds in this regime.

For pulse durations larger than 0.1 ms the LIDT for the irradiance has to also be taken into consideration!

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Pulse repetition rate R

Pulsed lasers with a high repetition rate may show very similar behavior to beams of cw lasers. As this is highly dependent on absorption and heat dissipation, there is no solid rule to determine when an optic will be damaged due to thermal effects. If you are unsure about the power level, please contact SILL.

Laser spot size (diameter d, or area A)

The spot size has an enormous impact on the resistance of the lens materials. As an example, assume a 1064 nm laser with a pulse energy of 10 mJ. When used with 500 μ m beam diameter through the lens, it would have a fluence of F=5.093 J/cm² and therefore surely destroy any material with the stated LIDT = 1 J/cm² from above. If the beam is expanded to a 5 mm beam diameter instead, the fluence would be F=0.051 J/cm² and the material should not be damaged at all.

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Pulse shape

The rule of thumb about the pulse shape is as follows: As peak intensities are considered, there is a difference between a top-hat like shape and a Gaussian beam. Just note that the LIDT value for a Gaussian beam is half the value of a top-hat shaped one.

Quick reference formula

For a quick estimate of the LIDT for your system with a laser pulse duration larger than 30 ps, the following formula can be used, while "spec" indicates the values from the specification:

E[J] or	P[W] > I IDT	$\tau \sim \frac{\lambda}{\tau}$, $\frac{\tau}{\tau}$, IDT	
$\left(\frac{d_{\rm L}}{2}\right)^2 \cdot \dot{\rm A}$	$\left(\frac{d_{\rm L}}{2}\right)^2 \cdot \dot{\rm A}$	$\sim \lambda_{\rm spec} \sqrt{\tau_{\rm spec}}$	

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As an example, a pulse laser at 532 nm wavelength and 50 ns pulse duration is used. The given LIDT from the lens supplier may again be 1 J/cm², measured at 1064 nm with a pulse duration of 10 ns. As calculated below, the LIDT can be estimated to be nearly 1.12 times higher.

 $\frac{532\,\text{nm}}{1064\,\text{nm}} \cdot \sqrt{\frac{50.0\,\text{ns}}{10.0\,\text{ns}}} \cdot 1.00 \frac{\text{J}}{\text{cm}^2} \approx 1.12 \frac{\text{J}}{\text{cm}^2}$

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Destruction mechanisms

The main damage mechanisms strongly depend on the laser properties, but are also dependent on the lens and coating materials in respect to geometrical properties, thermal conduction, absorption, defects, homogeneity, and more.

Thermal absorption has to be considered with the use of cw or long pulse (μ s) lasers. The damage is based on heating the lens material until melting, irreversible material expansions or modifications to the lens material. This damage mechanism also has to be considered for short-pulse lasers with high repetition rates, when the material is not given enough time to cool down between pulses.

Another cause of damage initially comes to play with short pulse lasers. Local defects or even tiny inclusions can cause plasma build up due to poor thermal conduction. This results in stress based breakthroughs and irreversible changes of the lens optical properties.

Furthermore, when considering ultra-short pulse lasers, additional multiphoton processes can occur ionizing coating materials and producing electrical breakthroughs when reaching critical electron densities.

LIDT of our coatings

In the following table we present the results of measurements of the laser induced damage threshold. They were done in a clean lab environment with acetone cleaned samples. Please be aware that due to coating process variations, geometry of the lens elements and the special testing environment, the real LIDT values might be much lower (rule of thumb is 5 to 10 times to be sure). Therefore these values are only test results and not a specification. In case of doubt, please contact our experienced team. Otherwise we cannot recompense in the event of damage.

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Correct lens handling and cleaning

To ensure a good performance and a long lifetime of optical components, the proper handling and cleaning of lens is critical. Contaminations such as dust, water and processing remnants can increase scatter light and light absorption and lead to defects. Thus, the optics should be stored in the delivered package and only opened in a clean environment. To avoid fingerprints on the optical surface and corresponding lens damage through skin oils, gloves should be worn. For a better handling of small lenses, the lens mount can be held with a tweezer. Take special care to avoid scratches on the lens surface.

If contamination is visible on the front lens of an optical system or on a single lens, it can be removed by these simple cleaning methods Dust can be removed by using inert dusting gas or a blower bulb. If there is still dust left on the surface, please use special lens tissues for cleaning. As required, the lens tissues can be soaked with acetone for a high cleansing effect and prevention of striations. With severe contamination, it may be necessary to clean the surface first with distilled water and if necessary, acetone.

Independent demounting of lenses leads to the expiration of the warranty and should be done only by the manufacturer!

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